

Are you considering investing in new benchtop equipment?

CORNING

Order Corning®, Axygen®, Benchtop Laboratory Equipment and Get a Special discount.

Corning® Lambda™ Plus Pipettor

- » Factory calibrated to strict ISO 8655 standards
- » Fully autoclavable and UV resistant
- » Designed for universal fit with all common brands of pipet tips
- » Adjustable tip ejector



Corning Lambda Plus

Cat. No.	Description	Qty/ Pk	Qty/ Cs	List Price	-%	Promo Price
4070	Corning Lambda Plus Single Channel Pipettor, 0.2 to 2 µL, Autoclavable	1	1	171,70 €	-20%	137 €
4071	Corning Lambda Plus Single Channel Pipettor, 0.5 to 10 µL, Autoclavable	1	1	171,70 €	-20%	137 €
4072	Corning Lambda Plus Single Channel Pipettor, 2 to 20 µL, Autoclavable	1	1	171,70 €	-20%	137 €
4073	Corning Lambda Plus Single Channel Pipettor, 10 to 100 µL, Autoclavable	1	1	171,70 €	-20%	137 €
4074	Corning Lambda Plus Single Channel Pipettor, 20 to 200 µL, Autoclavable	1	1	171,70 €	-20%	137 €
4075	Corning Lambda Plus Single Channel Pipettor, 100 to 1000 µL, Autoclavable	1	1	171,70 €	-20%	137 €
4080	Corning Lambda Plus 8 Channel Pipettor, 1 to 10 µL, Autoclavable	1	1	444,40 €	-25%	329 €
4081	Corning Lambda Plus 8 Channel Pipettor, 5 to 50 µL, Autoclavable	1	1	444,40 €	-25%	329 €
4082	Corning Lambda Plus 8 Channel Pipettor, 20 to 200 µL, Autoclavable	1	1	444,40 €	-25%	329 €
4083	Corning Lambda Plus 8 Channel Pipettor, 50 to 300 µL, Autoclavable	1	1	444,40 €	-25%	329 €
4084	Corning Lambda Plus 12 Channel Pipettor, 1 to 10 µL, Autoclavable	1	1	474,70 €	-25%	356 €
4085	Corning Lambda Plus 12 Channel Pipettor, 5 to 50 µL, Autoclavable	1	1	474,70 €	-25%	356 €
4086	Corning Lambda Plus 12 Channel Pipettor, 20 to 200 µL, Autoclavable	1	1	474,70 €	-25%	356 €
4087	Corning Lambda Plus 12 Channel Pipettor, 50 to 300 µL, Autoclavable	1	1	474,70 €	-25%	356 €
4069	Corning Lambda Plus Starter Kit	1 kit	1 kit	690,10 €	-25%	506 €

Corning® Stripettor™ Ultra Pipet Controller

- » Ergonomic, comfortable design
- » Variable speed range and pressure sensitive control buttons
- » Settings for blow out or slow « gravity » dispensing
- » Compatible with pipets from 1 to 100 mL



Corning Stripettor Ultra

Cat. No.	Description	Qty/Pk	Qty/Cs	List Price	-%	Promo Price
4099	Corning Stripettor Ultra, with Universal Power Supply, with 3 Filters 0,2 µm	1	1	290 €	- 30%	200 €

Corning® LSE™ Mini Microcentrifuge

- » Quick spin downs of micro-samples at a fixed speed of 6,000 rpm (2,000 x g)
- » Operation is simple and convenient
- » Compact & lightweight
- » Includes 8 x 1.5/2.0 mL rotor & 4-place strip tubes rotor



Corning LSE Equipment

Cat. No.	Description	Qty/Pk	Qty/Cs	List Price	-%	Promo Price
6770	Corning LSE Mini Microcentrifuge, 110-240V	1	1	195 €	-30%	140 €

Corning® LSE™ Compact Centrifuge

- » Capacity: 12 x 15 mL or 6 x 50 mL rotor
- » Accommodates also a swing-out rotor for 6 x 5 mL
- » Maximum speed: 6,000 rpm/4,185 x g
- » Accepts both conical and round bottom tubes



Corning LSE Equipment

Cat. No.	Description	Qty/Pk	Qty/Cs	List Price	-%	Promo Price
6759	Corning LSE Compact Centrifuge, 230 V EU PLUG	1	1	1650 €	-25%	1215 €

Axygen® Mini Plate Spinner Centrifuge

- » Maximum capacity: 2 standard PCR or microtiter plates
- » Ambient operating range 4° to 35°C
- » Top loading design
- » Accepts skirted and non-skirted and all standard PCR plates



Axygen Mini Plate Spinner Centrifuge

Cat. No.	Description	Qty/Pk	Qty/Cs	List Price	-%	Promo Price
Platespinner-230EU	Axygen® Mini Plate Spinner Centrifuge, 230V, EU Plug	1	1	595 €	-40%	350 €

Corning® LSE™ Vortex Mixer with Standard Tube Head

- Variable speed for gentle to vigorous mixing
- Optimized counter balance for minimal vibration
- Lightweight
- Includes cup head for all standard microtubes and centrifuge tubes, optional accessory heads available



Corning LSE Equipment

Cat. No.	Description	Qty/Pk	Qty/Cs	List Price	-%	Promo Price
6776	Vortex mixer with standard tube head, 230-V EU Plug	1	1	196,95€	-25%	147 €

Corning® Stirring Hot Plate

- Pyroceram® top
- Both temperature and stirring speed displayed on separate digital displays
- Temperature accuracy within $\pm 2^{\circ}\text{C}$ and stirring speed accuracy within $\pm 5\%$ of the speed setting chosen
- Digital display indicates when set temperature and stirring speed are reached



Corning Stirring Hot Plate

Cat. No.	Description	Qty/Pk	Qty/Cs	List Price	-%	Promo Price
6796-420D	Corning 5 x 7 Inch Top PC-420D Stirring Hot Plate/Stirrer with Digital Displays	1	1	535 €	-25%	395 €

Corning® LSE™ Digital Water Bath

- Compact design
- Stainless steel chamber and cover
- Powerful 500W heater programmed by LED digital controller ▸ Precise temperature control by hydraulic thermostat



Corning LSE Equipment

Cat. No.	Description	Qty/Pk	Qty/Cs	List Price	-%	Promo Price
6784	Corning LSE Digital Water Bath, 6L, 230V EU&UK Plug	1	1	656,50 €	-24%	500 €

Corning® LSE™ Digital Dry Bath Heater

- » Microprocessor control with large digital display
- » Temperature accuracy of $\pm 0.3^{\circ}\text{C}$ and uniformity of $\pm 0.2^{\circ}\text{C}$
- » Wide temperature set range from ambient $+5^{\circ}\text{C}$ to 150°C
- » Corrosion-resistant block cavity



Corning LSE Equipment

Cat. No.	Description	Qty/Pk	Qty/Cs	List Price	-%	Promo Price
6876-SB	Digital Dry Bath Heater, single block, 230V-EU plug	1	1	358,55 €	-24%	271 €
6886-DB	Digital Dry Bath Heater, dual block, 230V-EU plug	1	1	424,20 €	-25%	315 €
6896-FB	Digital Dry Bath Heater, four block, 230V-EU plug	1	1	626,20 €	-25%	471 €

Corning® LSE™ Shaking Incubators

- » Suitable for cold room or incubator use
- » Single or double knob control
- » Microplate shaker –process up to 4 microplates
- » Low Speed Orbital shaker- 19 mm orbit for blotting, gel staining and general mixing
- » Orbital shaker- 19 mm orbit for aeration of cultures



Corning LSE Equipment

Cat. No.	Description	Qty/Pk	Qty/Cs	List Price	-%	Promo Price
6791	Corning LSE Benchtop Shaking Incubator W/Platform, 230V EU&UK Plug	1	1	2990 €	-18%	2430 €
6752	Corning LSE Shaking Incubator, 49 Liters, 230V (EU and UK Plug)	1	1	4140 €	-24%	3125 €
6754	Corning LSE Shaking Incubator, 71 Liters, 230V (EU and UK Plug)	1	1	4715 €	-25%	3544 €

Corning® Cell Counter

- » Fast – thanks to its online image processing
- » Accurate – thanks to its cloud-based machine learning algorithm
- » Low-cost – works with common reusable glass hemocytometer. No consumables required.



Corning Cell Counter

Cat. No.	Description	Qty/Pk	Qty/Cs	List Price
6749	Corning Cell Counter	1	1	3355 €

CoolCell® Cell Freezing Containers

- » Alcohol-free cell-freezing containers ensure standardized controlled-rate -1°C/minute cell freezing in a -80°C freezer
- » Proven for use with a variety of cell types including stem cells, primary cells and cell lines
- » CoolCell technology utilizes a thermo-conductive alloy and highly-insulative outer materials to control the rate of heat removal and make cell cryopreservation reproducible




CoolCell Freezing Container

Cat. No.	Description	Qty/Pk	Qty/Cs	List Price
432138	CoolCell® LX-4 pack, 4 colours, Cell Freezing Container, for 12 x 1 ml or 2 ml cryogenic vials	1	1	560,45 €
432001	CoolCell® LX, purple Cell Freezing Container, for 12 x 1 ml or 2 ml cryogenic vials	1	1	165,20 €
432002	CoolCell® LX, green Cell Freezing Container, for 12 x 1 ml or 2 ml cryogenic vials	1	1	165,20 €
432003	CoolCell® LX, orange Cell Freezing Container, for 12 x 1 ml or 2 ml cryogenic vials	1	1	165,20 €
432004	CoolCell® LX, pink Cell Freezing Container, for 12 x 1 ml or 2 ml cryogenic vials	1	1	165,20 €
432005	CoolCell® 5 ml LX, purple Cell Freezing Container, for 12 x 3.5 ml or 5 ml cryogenic vials	1	1	221,50 €
432006	CoolCell® FTS30, purple Freezing Container, for 30 x 1 ml or 2 ml cryogenic vials	1	1	451,60 €
432007	CoolCell® FTS30, orange Freezing Container, for 30 x 1 ml or 2 ml cryogenic vials	1	1	451,60 €
432008	CoolCell® FTS30, green Freezing Container, for 30 x 1 ml or 2 ml cryogenic vials	1	1	451,60 €
432009	CoolCell® FTS30, pink Freezing Container, for 30 x 1 ml or 2 ml cryogenic vials	1	1	451,60 €

Axygen® Gel Documentation Systems

- » Auto exposure—you never have to manually manipulate the camera
- » Ready-to-use—just connect the computer and load the software
- » Long-life dual wavelength transilluminators
- » Exceptionally high resolution (5.4 MP) camera
- » The cabinet is compact with an ultra-violet (UV) cut-off switch to prevent exposure to UV light if the cabinet door is opened while the UV light is on



Axygen Gel Documentation systems

Cat. No.	Description	Qty/Pk	Qty/Cs	List Price	-%	Promo Price
GD-1000	Axygen Gel Documentation System	1	1	4950 €	-10%	4455 €
GDBL-1000	Axygen Gel Documentation System BL	1	1	8110 €	-7%	7513 €
HGB-7	Axygen horizontal gel box, 7 cm	1	1	309 €	-27%	225 €
HGB-10	Axygen horizontal gel box, 10 cm	1	1	365,65 €	-26%	269 €
HGB-15	Axygen horizontal gel box, 15 cm	1	1	551,05 €	-27%	400 €
HGB-20	Axygen horizontal gel box, 20 cm	1	1	618 €	-31%	425 €

Axygen® MaxyGene™ II Thermal Cycler

- » Unique flexible programming
- » Rapid run times
- » Improved workflow over traditional gradient cyclers
- » Ramping rates up to 5°C/sec
- » Simple user interface
- » Adjustable heated lid accommodates strips, tubes and microplates
- » Pre-programmed standard protocols



Axygen MaxyGene II Thermal Cycler

Cat. No.	Description	Qt/ Pk	Qty/ Cs	List Price	-%	Promo Price
THERM-1000	Axygen MaxyGene II Thermal Cycler with 96 well block, 230V	1	1	3846 €	-12%	3374 €

Axygen® PlateMax Semi Automated Plate Sealer

- » Adjustable sealing temperature: 100°C to 190°C
- » Compatible plate materials: PP (Polypropylene), PS (Polystyrene), PE (Polyethylene)
- » Compatible plate types: standard assay plates, deep well storage plates, PCR plates: skirted, semi-skirted, non-skirted



Axygen PlateMax

Cat. No.	Description	Qty/Pk	Qty/Cs	List Price	-%	Promo Price
HS-1230	Axygen PlateMax Semi Automated Plate Sealer, includes adapters for standard and deepwell microplates, 230V	1	1	4646.65 €	-35%	3000 €



Offer valid for purchases from August 1, 2018 through December 31, 2018.

Send your order by email: piia@inbio.ee

- » Put promotional code **EQUIP** on your purchase order
- » Apply any of the promotional prices above to the selected items
- » All prices marked here are without **VAT**. Transport cost will be not added.

Account Manager: Piia Kivipõld

Email: piia@inbio.ee

Mobile: +372 5301 6661

Offer valid in Estonia, Latvia and Lithuania only and void where prohibited by government, law, or company policy. Offer valid for end-user customers only for purchases from August 1, 2018 through December 31, 2018, or while supply lasts. For more specific information on claims, visit the Certificates page at www.corning.com/lifesciences. **Warranty/Disclaimer:** Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

Corning Incorporated

CORNING | **FALCON**® | **AXYGEN**® | **GOSSELIN**®

Corning BV
Fogstraat 12
1060 LJ Amsterdam
The Netherlands
Phone: +31 (0) 20 659 60 51
CSEurope@corning.com
corning.com/lifesciences
www.inbio.ee
Life Sciences Europe

For additional product or technical information, visit www.corning.com/lifesciences, or contact our Scientific Support Team at ScientificSupportEMEA@corning.com.

For a listing of trademarks, visit www.corning.com/clstrademarks. All other trademarks are the property of their respective owners.